

Department of Biology, Chemistry,  
Pharmacy  
Institute for Pharmacy  
(Pharmacology and Toxicology)  
Königin-Luise-Str. 2+4  
14195 Berlin  
Germany

## Standard Operating Procedure (SOP)

<b>Title</b> Viability assessment of cell monolayers with MTT reduction assay in 96-well plates			<b>Date</b> 2017-01-05	
<b>Document No.</b> SOP_MTT96		<b>First edition</b> 2014-07-03	<b>Issued by</b> N. Zhang, V. Kral	
		<b>Version 4</b>		
<b>Version</b>	<b>valid from</b>	<b>Description of changes</b>		
1	2014-07-03	Creation		
2	2014-12-01	Adjustment of serum concentration		
3	2015-10-15	Supplementary characterizations of nanocarriers; Adjustment of positive controls		
4	2017-01-05	Translation		
<b>Issued by</b> N. Zhang, V. Kral		<b>Reviewed by</b> C. Zoschke, C. Gerecke		<b>Approved by</b> Schäfer-Korting, M.
<b>Scope</b> Workgroups of the CRC112-Z01 and BB3R projects				

---

**Directory**

<b>1</b>	<b>Aims</b> .....	<b>3</b>
<b>2</b>	<b>Scope</b> .....	<b>3</b>
<b>3</b>	<b>Introduction</b> .....	<b>3</b>
<b>4</b>	<b>Material</b> .....	<b>3</b>
	Equipment .....	3
	Cells .....	5
	Consumables .....	6
	Cell Culture Medium .....	7
	Solutions .....	8
<b>5</b>	<b>Step-to-step protocol</b> .....	<b>9</b>
	Seeding Cells .....	9
	Nanocarrier Stimulation .....	9
	MTT Evaluation .....	10
<b>6</b>	<b>Data Analysis</b> .....	<b>11</b>
<b>7</b>	<b>Accepting Criteria</b> .....	<b>12</b>
<b>8</b>	<b>Annex</b> .....	<b>13</b>

<b>Issued by</b> V. Kral, N. Zhang	<b>Reviewed by</b> C. Gerecke, C. Zoschke.	<b>Approved by</b> Schäfer-Korting, M.	<b>Document No.</b> SOP_MTT96	<b>Date</b> 2017-01-05	<b>Vers.</b> 4
---------------------------------------	---	---	----------------------------------	---------------------------	-------------------

## AIMS

Standardized viability assessment of cell monolayers with MTT (3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide) reduction assay for the cytotoxic potential of test nanocarriers.

## SCOPE

This SOP applies to the workgroups of CRC112 and BB3R projects.

## INTRODUCTION

The MTT assay is applied to assess the viability of cell monolayers after treated with test nanocarriers. The method is based on the colorimetric change of the yellow tetrazolium dye into the purple formazan by oxidoreductase enzymes in viable cells.

For the study of nanocarriers, the MTT assay is performed on normal human keratinocytes (NHK) and normal human dermal fibroblasts (NHDF) isolated from foreskins. The primary cells are pooled after isolation (3 donors per cell type).

The physicochemical characterization of the nanocarriers is provided in data sheets, see Annex 1, by project partners, such as Prof. Bodmeier, Prof. Haag, Prof. Calderon and Prof. Lendlein, etc.

## MATERIALS

### Equipment

Designation	Manufacturer
Autoclave	Systemec, Wetzlar
Cell counting chamber	Zeiss, Jena
Centrifuge (Eppendorf)	Eppendorf, Hamburg
Centrifuge (Megafuge <sup>®</sup> 1.0R)	Thermo Fisher Scientific, Waltham, MA, USA
Fluorescence microscope (BZ-8000K)	Keyence, Osaka, JAP
Fluostar Optima	BMG Labtech, Offenbach
Freezer (-20°C)	Siemens, München
Freezer (-80°C)	Thermo Fisher Scientific, Waltham, MA, USA
Incubator (BB6220)	Thermo Fisher Scientific, Waltham, MA, USA
Magnetic stirrer RCT basic	IKA-Werke, Staufen
Microtome (Hyrax M40)	Zeiss, Jena
Nitrogen tank (Arpege 70)	Air Liquide, Paris, F
Phase contrast inverted microscope (Axiovert 40C)	Zeiss, Jena

Issued by	Reviewed by	Approved by	Document No.	Date	Vers.
V. Kral, N. Zhang	C. Gerecke, C. Zoschke	Schäfer-Korting, M.	SOP_MTT96	2017-01-05	4

---

pH meter (766 Calimatic)	Knick, Nürnberg
Pipette (Eppendorf Reference®)	Eppendorf, Hamburg
Pipetting aid (Easypet®)	Eppendorf, Hamburg
Refrigerator (4°C)	Siemens, München
Shaker IKA® MTS 2	IKA, Staufen
Sterile working bench (LaminAir®)	Thermo Fisher Scientific, Waltham, MA, USA
Suction pump	VWR, Darmstadt
Vortex	Bender & Hobei, Zurich, CH
Water bath	Gesellschaft für Labortechnik, Burgwedel
Water processing unit (SG LaboStar)	SG Wasseraufberitung und Regenerierstation, Barsbüttel

---

Equivalent equipment of other suppliers can be used as well.

<b>Issued by</b> V. Kral, N. Zhang	<b>Reviewed by</b> C. Gerecke, C. Zoschke.	<b>Approved by</b> Schäfer-Korting, M.	<b>Document No.</b> SOP_MTT96	<b>Date</b> 2017-01-05	<b>Vers.</b> 4
---------------------------------------	---	---	----------------------------------	---------------------------	-------------------

---

**Cells**

<b>Designation</b>	<b>Source</b>
NHDF, passage 3, pooled of 3 donors	Isolation from juvenile preputium <sup>1)</sup>
NHK, passage 3, pooled of 3 donors	Isolation from juvenile preputium <sup>1)</sup>

**Abbreviations**

<b>Designation</b>	<b>Supplier</b>
5-FU	5-fluorouracil
BPE	Bovine pituitary extract
DMEM	Dulbecco's modified Eagle's medium
DMSO	Dimethyl sulfoxide
EDTA	Ethylenediaminetetraacetic acid
FCS	Fetal calf serum
KCl	Potassium chloride
KH <sub>2</sub> PO <sub>4</sub>	Potassium dihydrogenphosphate
MTT	3- (4,5-Dimethylthiazol-2-yl) -2,5-diphenyltetrazolium bromide
NaCl	Sodium chloride
Na <sub>2</sub> HPO <sub>4</sub>	Disodium hydrogen phosphate
OD	Optical density
PBS	Phosphate buffered saline solution (pH 7.4)
SDS	Sodium dodecyl sulfate

---

<sup>1)</sup> according to SOP "Isolation of keratinocytes and fibroblasts from human specimens"

<b>Issued by</b>	<b>Reviewed by</b>	<b>Approved by</b>	<b>Document No.</b>	<b>Date</b>	<b>Vers.</b>
V. Kral, N. Zhang	C. Gerecke, C. Zoschke.	Schäfer-Korting, M.	SOP_MTT96	2017-01-05	4

## Consumables

Designation	Supplier
96-well plate (flat bottom)	TPP, Trasadingen, Switzerland
Cell culture flask (75 cm <sup>2</sup> )	TPP, Trasadingen, Switzerland
Cell culture flask (150 cm <sup>2</sup> )	TPP, Trasadingen, Switzerland
Centrifuge tube (15 mL)	Sarstedt, Nümbrecht
Centrifuge tube (50 mL)	Sarstedt, Nümbrecht
DMEM	Sigma-Aldrich, München
DMSO	Sigma-Aldrich, München
EDTA	Sigma-Aldrich, München
FBS Superior	Biochrom, Berlin
H <sub>2</sub> O ( pyrogen free)	Carl Roth, Karlsruhe
KCl	Sigma-Aldrich, München
KH <sub>2</sub> PO <sub>4</sub>	Carl Roth, Karlsruhe
KGM Bulletkit	Lonza, Köln
L-Glutamin	Biochrom, Berlin
MTT	Sigma-Aldrich, München
Na <sub>2</sub> HPO <sub>4</sub>	Carl Roth, Karlsruhe
NaCl	Carl Roth, Karlsruhe
Keratinocytes Growth Medium (KGM)	Lonza, Köln
L-Glutamine	Sigma-Aldrich, München
Penicillin-Streptomycin-solution	Biochrom, Berlin
Pipette tips	Eppendorf, Hamburg
Trypan blue	Biochrom, Berlin
Trypsin dry substance	Biochrom, Berlin

Equivalent consumables of other suppliers can be used as well.

Issued by	Reviewed by	Approved by	Document No.	Date	Vers.
V. Kral, N. Zhang	C. Gerecke, C. Zoschke.	Schäfer-Korting, M.	SOP_MTT96	2017-01-05	4

### Cell culture medium

Description	Ingredients	Remarks
FBM	DMEM	500 mL
	L-Glutamine	5 mL
	Penicillin/Streptomycin	5 mL
FGM	FBM	510 mL
	Fetal Calf Serum	37.5 mL
KGM	KBM	500 mL
	KGM supplements	

For NHDF, FGM is applied as cell culture medium, and FBM as nanocarrier testing medium.

For NHK, KGM is applied as both cell culture and nanocarrier testing medium.

### General information

Unless stated otherwise, all procedures should be performed under sterile laminar flow conditions, cells are incubated in the incubator at 37°C, 5 % CO<sub>2</sub> and 95 % humidity, and all the mediums are pre warmed in water bath at 37°C for 30 min.

### Colored nanocarriers

For assessment of colored nanocarriers, before undergoing the procedure, first mix the testing nanocarrier at 0.05% (w/v) with freshly prepared MTT solution. If the solution turns blue or purple, the nanocarriers are supposed to reduce MTT directly.

And to further study the interference of the colorant with the MTT reading, an additional color control should be performed. When performing the standard absorbance measurement, each interfering nanocarrier is applied in triplicates per exposure time, which undergoes the same procedures but is incubated with medium instead of MTT solution during the MTT incubation step to generate a non-specific color control. The final viability of the test nanocarriers is then calculated by subtracting the value of the color control.

When the MTT reading does not meet the accepting criteria, an alternative method, e.g. Neutral Red Uptake assay (OECD Test Guideline 432), could be performed for the aim of cytotoxicity evaluation.

Issued by	Reviewed by	Approved by	Document No.	Date	Vers.
V. Kral, N. Zhang	C. Gerecke, C. Zoschke.	Schäfer-Korting, M.	SOP_MTT96	2017-01-05	4

## Solutions

### PBS and Trypsin-EDTA solutions

PBS and Trypsin-EDTA solutions are prepared according to the SOP as on page 5, autoclaved, and stored at 4°C.

### Positive Control SDS solutions

Dilute 5 mg of SDS with 10 mL of distilled H<sub>2</sub>O into a 0.05% (w/v) stock solution, and a 1:10 dilution of the stock solution and respective testing medium is applied as positive control.

### Solvent Control solutions

Dilute pyrogen free H<sub>2</sub>O or PBS, depending on the solvent used in the nanocarrier respectively, in 1:10 with respective testing medium. The solution is applied as solvent control,.

### MTT solutions

Dilute MTT powder with PBS into a 5 mg/mL solution, aliquot and store at -20°C as stock solution. Dilute MTT stock solution 1:10 with respective testing medium.

### Testing nanocarriers solutions

Dilute the nanocarrier solution with pyrogen free H<sub>2</sub>O or PBS, depending on the solvent used in the nanocarrier respectively, into 0.5% and 0.05% (w/v) solutions. Dilute these nanocarriers 1:10 with the respective testing medium to a final concentrations of 0.05% and 0.005% (w/v).

<b>Issued by</b> V. Kral, N. Zhang	<b>Reviewed by</b> C. Gerecke, C. Zoschke.	<b>Approved by</b> Schäfer-Korting, M.	<b>Document No.</b> SOP_MTT96	<b>Date</b> 2017-01-05	<b>Vers.</b> 4
---------------------------------------	---	---	----------------------------------	---------------------------	-------------------



## STEP-TO-STEP PROTOCOL

### Seeding cells | Day 0

- 1.5 mL of Trypsin-EDTA is added into the cell culture flask 75 cm<sup>2</sup> and put into the incubator for 3 min.
- 3.5 mL of KGM or FGM is added into the flask, and the cell suspension is transferred into a 50 mL centrifuge tube. The cell culture bottom is then washed twice with 5 mL PBS each, and all washed cell suspensions are transferred to the same tube
- Centrifugate at 1000 rpm for 5 min. discard the supernatant and resuspend the cells with 10 mL of PBS, then count the cells according to the SOP on page 5 and repeat the centrifugation.
- Adjust the cell suspension to  $1 \times 10^5$  cells / mL with relative cell culture medium.
- Pipet 100  $\mu$ L of the cell suspension into each inner wells of the 96-well plate (see scheme day 0), which corresponds to 10,000 cells / well. And pipet 100  $\mu$ L of the relative cell culture medium, see Page 7, into the outer wells. The plates are then placed into the incubator for 24 h.

	1	2	3	4	5	6	7	8	9	10	11	12
A	B	B	B	B	B	B	B	B	B	B	B	B
B	B	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	B
C	B	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	B
D	B	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	B
E	B	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	B
F	B	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	B
G	B	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	B
H	B	B	B	B	B	B	B	B	B	B	B	B

#### **Scheme day 0: cell seeding**

*B: Blank (Cell culture medium); CS: cell suspension ( $1 \times 10^5$  cells / ml)*

### Nanocarrier Stimulation | Day 1

- Prepare all the solutions as described in the solutions section.
- Check the cell morphology under the microscope.
- Aspirate all the medium from the plate with the suction pump and pipet tips (yellow) in a gentle manner.

<b>Issued by</b> V. Kral, N. Zhang	<b>Reviewed by</b> C. Gerecke, C. Zoschke.	<b>Approved by</b> Schäfer-Korting, M.	<b>Document No.</b> SOP_MTT96	<b>Date</b> 2017-01-05	<b>Vers.</b> 4
---------------------------------------	---	---	----------------------------------	---------------------------	-------------------

- 100  $\mu$ L of the corresponding medium or solutions are then pipetted into each well (see scheme day 1).
- The plates are then placed into the incubator for 24 h or 48 h.

	1	2	3	4	5	6	7	8	9	10	11	12
A	B	B	B	B	B	B	B	B	B	B	B	B
B	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
C	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
D	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
E	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
F	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
G	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
H	B	B	B	B	B	B	B	B	B	B	B	B

**Scheme day 1**

*B: Blank (only testing medium without cells); UC: Untreated Control (only testing medium with cells); SC: Solvent Control; PC: Positive Control; C1-6: Testing nanocarrier solutions (a: 0.05%, b: 0.005%)*

### MTT Evaluation | Day 02 or 3

- Prepare the MTT solution as described in the solutions section.
- Check the cell morphology under the microscope.
- Aspirate all the medium from the plate with the suction pump and pipet tips (yellow) in a gentle manner, and wash each well once with 100  $\mu$ L of PBS.
- Add 100  $\mu$ L of the MTT solution into each well, then place the plates into the incubator for 4 h.
- Remove the MTT solution by gentle suction and place the plates over a sterilized paper towel for 1 min.
- Add 50  $\mu$ L of DMSO into each well, and put the plate onto the plate shaker at 500 rpm for 10 min.
- Measure the absorption at the wavelength of 540 nm, name all the relative excel and prism files with the same experiment ID, and calculate the viability.

### Photometer – Setting

Device	FLUOstar OPTIMA (or equivalent device)
Programme	MTT Absorbance, TPP96
Mode	Disk Mode
Positioning Delay	0.7 s
No. Kinetic Windows	1

Issued by	Reviewed by	Approved by	Document No.	Date	Vers.
V. Kral, N. Zhang	C. Gerecke, C. Zoschke.	Schäfer-Korting, M.	SOP_MTT96	2017-01-05	4

---

Excitation Filter	A540
Emission Filter	Empty
No. Multichromatics	1
Shaking width	1 mm; 600 rpm
No. Cycles	1
No. Flashes per well and cycle	2

---

A comparable device with adequate adjustment can also be used.

### Data Analysis

- The blank value is subtracted from all measured values to gain a corrected OD value.
- From each plate, the mean corrected value for the solvent control is set equal to 100%. The viability rates of the test nanocarriers are calculated as follows:

$$\text{Viability (\%)} = \frac{\text{Corrected Testing Nanocarrier Value}}{\text{Corrected Solvent Control}} \times 100\%$$

- An Excel and GraphPad Prism 5 spreadsheet is provided which is used for the calculation of the results. For use details of the spreadsheet see Annex 2.

<b>Issued by</b> V. Kral, N. Zhang	<b>Reviewed by</b> C. Gerecke, C. Zoschke.	<b>Approved by</b> Schäfer-Korting, M.	<b>Document No.</b> SOP_MTT96	<b>Date</b> 2017-01-05	<b>Vers.</b> 4
---------------------------------------	---	---	----------------------------------	---------------------------	-------------------

## Accepting Criteria

The results are acceptable if:

- The corrected OD values of the untreated controls are within the range:

Cell Type	Corrected OD Value
NHK	0.15–1.0
NHDF	0.3–1.0

- The difference of viability among triplicates is  $< 20\%$  in the same run.
- The values of the viability from column 2 and 11 (untreated controls) show a deviation of  $\leq 15\%$ .
- The viability of positive controls is below  $15\%$ .

Issued by	Reviewed by	Approved by	Document No.	Date	Vers.
V. Kral, N. Zhang	C. Gerecke, C. Zoschke.	Schäfer-Korting, M.	SOP_MTT96	2017-01-05	4

## Annex 1

Sample Data Sheet

Working Group \_\_\_\_\_

Project Name \_\_\_\_\_

## Table 1 General Information

<b>Sample Name</b>	<b>Date of Delivery</b>
<b>Technician or PhD Name</b>	<b>Tel. No.</b>
<b>Lab Journal No.</b>	<b>Batch No.</b>

Chemical Structure, Formulation/Solvent, Guest, Marker:

--

## Table 2 Delivered Amount, Concentration and Solubility

<b>Amount (mg or mL)</b>
<b>C(carrier)</b>
<b>Dissolved in</b>
<b>Sterilization</b>
<b>Storage Recommendation</b>
<b>History of Sample</b>

## Table 3 Characterization (optional)

<b>IR</b>	<b>DLS</b>
<b>NMR</b>	<b>REM</b>
<b>UV</b>	<b>TEM</b>

<b>Issued by</b> V. Kral, N. Zhang	<b>Reviewed by</b> C. Gerecke, C. Zoschke.	<b>Approved by</b> Schäfer-Korting, M.	<b>Document No.</b> SOP_MTT96	<b>Date</b> 2017-01-05	<b>Vers.</b> 4
---------------------------------------	---	---	----------------------------------	---------------------------	-------------------

## Annex 2

Table 1 (Excel) raw data OD

Experiment ID:

	1	2	3	4	5	6	7	8	9	10	11	12
A	B	B	B	B	B	B	B	B	B	B	B	B
B	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
C	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
D	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
E	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
F	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
G	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
H	B	B	B	B	B	B	B	B	B	B	B	B

Table 2 (Excel) corrected OD (raw data OD – blank OD)

Experiment ID:

	1	2	3	4	5	6	7	8	9	10	11	12
A	B	B	B	B	B	B	B	B	B	B	B	B
B	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
C	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
D	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
E	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
F	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
G	B	UC	SC	PC	C1b	C2b	C3b	C4b	C5b	C6b	UC	B
H	B	B	B	B	B	B	B	B	B	B	B	B
MEAN OD	B	UC	SC	PC	C1a	C2a	C3a	C4a	C5a	C6a	UC	B
	B (1+12)	UC (2+11)			C1b	C2b	C3b	C4b	C5b	C6b		B

Issued by V. Kral, N. Zhang	Reviewed by C. Gerecke, C. Zoschke.	Approved by Schäfer-Korting, M.	Document No. SOP_MTT96	Date 2017-01-05	Vers. 4
--------------------------------	--	------------------------------------	---------------------------	--------------------	------------

**Table 3 (Excel) viability (%)**

<b>Sample Name</b>	<b>Corrected OD</b>	<b>Viability (%)</b>
SC		100
PC		
C1a		
C1b		
C2a		
C2b		
C3a		
C3b		
C4a		
C4b		
C5a		
C5b		
C6a		
C6b		


<b>Issued by</b> V. Kral, N. Zhang	<b>Reviewed by</b> C. Gerecke, C. Zoschke.	<b>Approved by</b> Schäfer-Korting, M.	<b>Document No.</b> SOP_MTT96	<b>Date</b> 2017-01-05	<b>Vers.</b> 4
---------------------------------------	---	---	----------------------------------	---------------------------	-------------------

Table 4 (Prism)

Graph Collum, Plot Mean  $\pm$  SEM

Experiment ID:

	C1b	C1a	C2b	C2a	C3b	C3a	C4b	C4a	C5b	C5a	C6b	C6a	PC	SC
	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
1														
2														
3														



Issued by	Reviewed by	Approved by	Document No.	Date	Vers.
V. Kral, N. Zhang	C. Gerecke, C. Zoschke.	Schäfer-Korting, M.	SOP_MTT96	2017-01-05	4